



"EZ Mold" Diabetic Inserts

Save Time and Money when Forming Custom Diabetic Insoles Perfect for O & P and Pedorthotic Professionals

Our EZ Mold Inserts are designed for forming CUSTOM Diabetic Insoles. Use like a foam blocker, vacuum forming over a patient mold, but with improved product features AND less trimming, grinding and waste.

FEATURES:

Tri-Layer Construction - 3D Version 1/8" 20 Shore A ProCell top layer, 1/6" Poron cushioning middle layer, over a firm 40 Shore A ProCell EVA base.

Full Length - Total contact insole for use on the diabetic foot and at risk patients.

Built in Arch Support - Pre-formed arch support provides excellent longitudinal support. Makes forming quick and easy, eliminating the need for arch filling, grinding and sanding in most cases.

Metatarsal Raise - Can be adjusted when heated to increase met raise or pushed out to lower for sensitive feet.

Pre-formed Heel Cup, Arch Support - Pre-shaped with a deep heel cup, arch support, and medial/lateral rims to help save time and facilitate easier custom molding.

High Quality Materials - With cushioning Poron middle layer for long lasting comfort and shock absorption.

Heat Moldable - Can be formed at 250° Fahrenheit with a heat gun or placed in oven.



1-888-894-4171

www.shercoinc.com

Size	Women	Men
36	5 - 6	
38	7 - 8	
40	9 - 10	6 - 7
42	11 - 12	8 - 9
44	13	10 - 11
46		12 - 13
48		14+

Call Your Local Representative:

Discover the Entire Sher Brand Line of Professional Medical Orthotics and Materials.